

## Slide-Warmers

### Our Large Size Economical Slide-Warmer!

- Temperature range from slightly above ambient to 65°C
- Two sizes—accommodates up to 66 slides
- Large slide capacity with space-saving design
- ★ 24 Month Warranty

#### Product Description

- Uniform heat transfer over entire heating surface.
- Anodized black heavy gauge aluminum work top resists corrosion and provides high-contrast background for viewing specimen.
- Powder coated heavy-duty steel body construction is durable and easy to clean.
- Thermostat-controller regulates temperature to 65°C.
- Model 26025Q features LED display and cover.
- Model 26007CEQ and Model 26022CEQ are CE certified.
- Models 26005Q and 26020Q are ETL marked.
- Models 26005Q & 26007CE have a slide capacity of sixty-five 1" x 3" (2.54 cm x 7.6 cm) or thirty-two 2" x 3" (5.1 cm X 7.6 cm).
- Models 26020Q, 26022CEQ, 26025Q & 26025-1Q have a slide capacity of sixty-six 1" x 3" (2.54 x 7.6 cm) or thirty-three 2" x 3" (5.1 x 7.6 cm).

#### Application

- Designed for cytology, histology, and other biological applications in industrial and clinical lab environments.



Model 26025Q



Model 26005Q



#### PRODUCT SPECIFICATIONS

Model #	Electrical (50/60 Hz)		Overall Dimensions Inches (cm)		
	Volts	Watts	W	H	D
26005Q <sup>2</sup>	120	160	14 (36)	14 (36)	3.5 (9)
26007CEQ <sup>1</sup>	240	160	14 (36)	14 (36)	3.5 (9)
26020Q <sup>2</sup>	120	160	25 (64)	8 (20)	3.5 (9)
26022CEQ <sup>1</sup>	240	160	25 (64)	8 (20)	3.5 (9)
26025Q	120	160	25 (64)	8 (20)*	4.25 (11)
26025-1Q	240	160	25 (64)	8 (20)*	4.25 (11)

<sup>1</sup> CE

<sup>2</sup> ETL

\*Height with lid open

#### ACCESSORIES

Model #	Description
26023Q	Metal dust cover for models 26022 and 26020