


## Mini UV Viewing Cabinets

The mini Chromato-Vue® UV viewing cabinets are designed with a metal base and plastic top. The C-15G cabinet includes a UVGL-58 UV lamp. Lamps for other mini cabinets are ordered separately.

- The bottom panel of the cabinet is removable
- The lamp fits snugly on top of the cabinet
- A contrast control filter absorbs UV energy and protect the eyes
- Soft rubber viewport provides a comfortable viewing area with space for eyeglasses
- Lightweight access curtain blocks out external light from the cabinet interior



**C-10E Mini UV Cabinet**  
Designed for use with the 4-watt EL Series Lamps (see pages 21 - 22).



**C-15G Mini UV Cabinet**  
Combines the C-10 Cabinet with UVGL-58 multi-band shortwave/longwave lamp.

Cabinet Dimensions:

Interior: 9L x 10W x 7H in. (229 x 254 x 178mm)

Exterior: 9L x 10.5W x 12H in. (229 x 267 x 305mm)

### Ordering Information

<b>C-10</b> Mini Cabinet	95-0072-01
<b>C-10P</b> Portable Cabinet	95-0072-02
<b>C-10E</b> Mini Cabinet	95-0072-08
<b>C-15G</b> Cabinet with Lamp	95-0072-06 (115V) 95-0072-07 (230V)

**C-10P Mini UV Cabinet**  
accommodates the 6-watt Rechargeable Lamps which are excellent for field use applications (see page 24 for Rechargeable Lamps).



## Multifunctional UV Viewing Cabinet

C-65 Chromato-Vue Cabinet accommodates one or two EL Series 8-watt lamps ordered separately, (for additional lamps, see pp. 21, 22). Use of two lamps illuminates the cabinet's interior more uniformly from two sides of the cabinet.

- The bottom panel of the cabinet is removable
- Includes an interior overhead white light

### Ordering Information

<b>C-65</b> Mid-Size UV Viewing Cabinet	95-0257-01 (115V) 95-0257-02 (230V)
---	--

Dimensions C-65:

Int: 17L x 14W x 6H in. (533 x 330 x 400mm)

Ext: 19.75L x 16.5W x 15.75H in. (602 x 419 x 400mm)

EL Series eight-watt lamps are designed to be used with the C-65 cabinet or removed for handheld use.

- Lamps are made with a durable, scratch-resistant metal housing
- Lamps include a universal plug, which connects directly to the C-65 Cabinet, plus a supplementary power cord to plug into a line current



**C-65 UV Cabinet** is designed to accommodate one or two of the EL Series 8-watt Lamps

### Ordering Information - EL Lamps

<b>UVL-28</b> 365nm/365nm	95-0248-01 (115V) 95-0248-02 (230V)
<b>UVM-28</b> 302nm/302nm	95-0250-01 (115V) 95-0250-02 (230V)
<b>UVS-28</b> 254nm/254nm	95-0249-01 (115V) 95-0249-02 (230V)
<b>UVLM-28 (2UV)</b> 365nm/302nm	95-0251-01 (115V) 95-0251-02 (230V)
<b>UVLS-28 (2UV)</b> 365nm/254nm	95-0201-01 (115V) 95-0201-02 (230V)
<b>UVLMS-38 (3UV)</b> 365nm/302nm/254nm	95-0252-01 (115V) 95-0252-02 (230V)