

High Performance UV Viewing Cabinets



C-75 Large UV Cabinet features shortwave and longwave plus interior white light.

Dimensions C-75:

Interior: 17W x 14D x 6H in. (533 x 330 x 400mm)

Exterior: 19.75W x 16.5D x 15.75H in. (602 x 419 x 400mm)

Weight: 18 lb (8 kg)

C-70 and C-71 Cabinets

The C-70 and C-71 viewing cabinets feature a removeable bottom panel so the cabinet may be placed over large samples or a transilluminator. These cabinets are constructed of durable, lightweight plastic.

- The **C-70G** Chromato-Vue Cabinet includes two each of 15-watt shortwave (254nm) and longwave UV mounted on one side of the cabinet
- The **C-71** is similar to the C-70G, but has four 15-watt, longwave (365nm) UV tubes on one side of the unit
- Control switches are conveniently located for selecting either of the UV wavelengths or the interior overhead white light
- Easy access to the cabinet interior is via the curtain which also blocks ambient light
- The longwave provides brilliant fluorescence for applications such as non-destructive testing and document verification
- A soft rubber viewport includes a contrast control filter

Dimensions C-70G/C-71:

Interior: 21W x 13D x 8H in. (432 x 356 x 127mm)

Exterior: 23.4W x 15.8D x 13H in. (594 x 400 x 30mm)

Weight: 22 lb (10.3 kg)

C-75 Cabinet

High intensity 15-watt shortwave (254nm) and longwave (365nm) UV tubes are positioned on two sides of the C-75 Chromato-Vue Cabinet providing overhead ultraviolet illumination for many applications including fluorescence studies and non-destructive testing.

- Power switches are conveniently placed for operating the UV and interior white light
- A lightweight curtain blocks exterior light from the cabinet interior
- Soft rubber viewport provides a comfortable viewing area with adequate space for glasses

Certifications: All 115V models US/C Entela certified and 230V CE certified except where noted.



C-70 UV Cabinet includes two each of 15-watt short-wave and longwave UV tubes mounted on one side of the cabinet; the **C-71** has four 15-watt longwave tubes.

Ordering Information

| | |
|---|--|
| C-75 Large UV Viewing Cabinet 254/365nm UV | 95-0253-01 (115V) 95-0253-02 (230V) |
| C-70G Large Viewing Cabinet, 254/365nm UV | 95-0020-04 (115V) 95-0020-05 (230V) |
| C-71 Large Viewing Cabinet, 365nm UV | 95-0022-01 (115V) 95-0022-02 (230V) |

Replacement Parts

| | |
|---------------------|------------|
| 15-watt tube, 365nm | 34-0009-02 |
| 15-watt tube, 254nm | 34-0008-01 |