



## Now Available: Multi Shot 0.75kw System

5 cycles = 1 process (1-40mls per process)

Adjustable pressure – 1-40Kpsi

Reliable & continuous results – fixed orifice

Quick and easy to use

No need for compressed air or gas


Minimal dead loss (<1mls)

Heat sensitive samples are protected by cooling the collection cups prior to use

Suitable for most samples e.g. – bacteria, yeast, plant and mammalian tissue

### Can process the following sample states:

- Liquid
- Pellet
- Paste
- Solid
- Frozen



Footprint: 475mm x 510mm x 500mm (DxWxH)  
Weight: 110kgs

### High Pressure Cell Disruption

"A simple solution for cell disruption.  
Suitable for applications such as Nucleic  
acid extraction, protein extraction &  
tissue disruption."



### UK Headquarters

Constant Systems Limited,  
Low March, Daventry, Northants  
NN11 4SD  
ENGLAND  
UK

Tel: +44 (0) 1327 314146

Fax: +44 (0) 1327 314147

Email: [sales@constantsystems.com](mailto:sales@constantsystems.com)

[www.constantsystems.com](http://www.constantsystems.com)