

3UV™ Multi-Wavelength UV Lamps

This innovative lamp (patented) houses three ultraviolet wavelengths operated via A/C (mains) power.

- Selection of each wavelength is by an easy turn of the dial located on one end of the housing
- The uniquely designed reflectors mounted behind each tube provides maximum UV for fluorescence uses
- Lightweight and ergonomically designed for hand or stationary use

Ordering Information

3UV-34 4 Watts 254nm/302nm/365nm	95-0341-01 (115V) 95-0341-02 (230V)
3UV-36 6 Watts 254nm/302nm/365nm	95-0342-01 (115V) 95-0342-02 (230V)
3UV-38 8 Watts 254nm/302nm/365nm	95-0343-01 (115V) 95-0343-02 (230V)

For replacement tubes: see bottom of page.

EL Series UV Lamps

EL Series UV lamps are available in 4, 6, or 8-watt models and a selection of UV wavelengths and white light combinations. The lamps can be placed in the J-138 stand (see p. 22) for hands-free use. EL Series Lamps are manufactured with a scratch resistant powder painted metal housing.

The **3UV™** model is available with 8-watt tubes.

- 3UV model houses three UV wavelengths (shortwave 254nm, midrange 302nm, and longwave 365nm) in one unit reducing the need to purchase multiple lamps
- Easily switch among the three UV wavelengths as required for your application

The **2UV™** models include two different ultraviolet wavelengths in one lamp.

Lamp Dimensions:

4-Watt: 9.8L x 3.8W x 2.5D in. (249 x 97 x 64mm)
6-Watt: 12L x 3.8W x 2.5D in. (302 x 97 x 64mm)
8-Watt: 14.8L x 3.8W x 2.5D in. (376 x 97 x 64mm)


Replacement Tubes

Tube, 8-watt, White	34-0056-01
Tube, 8-watt, 365nm	34-0006-01
Tube, 8-watt, 302nm	34-0042-01
Tube, 8-watt, 254nm	34-0007-01
Tube, 6-watt, 365nm	34-0034-01
Tube, 6-watt, 302nm	34-0044-01
Tube, 6-watt, 254nm	34-0013-01
Tube, 4-watt, 365nm	34-0005-01
Tube, 4-watt, 254nm	34-0003-01

3UV Lamps
house longwave,
shortwave and
midrange UV in
one lamp
(8-watt model
pictured).

Lamp Dimensions:

3UV-34: 9.5L x 3W x 4.5D in. (241 x 76 x 114mm)
3UV-36: 12.5L x 3W x 4.5D in. (318 x 76 x 114mm)
3UV-38: 15.5L x 3W x 4.5D in. (394 x 76 x 114mm)



EL Series UV Lamps are lightweight and easy to handle. Left to right: UVLMS-38, UVLS-28, UVL-14.

Ordering Information

UVLMS-38 3UV 365nm/302nm/254nm Replacement Filter	95-0252-01 (115V) 95-0252-02 (230V) 38-0158-06
UVLM-28 8-watt 365nm/302nm 2UV Replacement Filter	95-0251-01 (115V) 95-0251-02 (230V) 38-0158-03
UVLS-28 8-watt 365nm/254nm 2UV Replacement Filter	95-0201-01 (115V) 95-0201-02 (230V) 38-0158-03
UVLM-26 6-watt 365nm/302nm 2UV Replacement Filter	95-0278-01 (115V) 95-0278-02 (230V) 38-0182-03
UVLS-26 6-watt 365nm/254nm 2UV Replacement Filter	95-0279-01 (115V) 95-0279-02 (230V) 38-0182-04
UVLS-24 4-watt 365nm/254nm 2UV Replacement Filter	95-0271-01 (115V) 95-0271-02 (230V) 38-0183-03