


EL Series UV Lamps

Basic and Twin EL Series Lamps are available in 4, 6 and 8 watt configurations.

- **Basic Models** include one UV tube (254nm, 302nm, or 365nm) and one white light tube in one lamp
- **Twin Models** include two tubes of the same UV wavelength for added UV intensity



EL Series Twin Models have two tubes of the same UV wave-length (UVL-28 shown).

Twin Models Ordering Information

UVL-28 8-watt, 365nm	95-0248-01 (115V) 95-0248-02 (230V)
Replacement Filter	38-0158-01
UVM-28 8-watt, 302nm	95-0250-01 (115V) 95-0250-02 (230V)
Replacement Filter	38-0158-01
UVS-28 8-watt, 254nm	95-0249-01 (115V) 95-0249-02 (230V)
Replacement Filter	38-0158-02
UVL-26 6-watt, 365nm	95-0275-01 (115V) 95-0275-02 (230V)
Replacement Filter	38-0182-05
UVM-26 6-watt, 302nm	95-0276-01 (115V) 95-0276-02 (230V)
Replacement Filter	38-0182-04
UVS-26 6-watt, 254nm	95-0277-01 (115V) 95-0277-02 (230V)
Replacement Filter	38-0182-04
UVL-24 4-watt, 365nm	95-0267-01 (115V) 95-0267-02 (230V)
Replacement Filter	38-0183-05
UVM-24 4-watt, 302nm	95-0268-01 (115V) 95-0268-02 (230V)
Replacement Filter	38-0183-04
UVS-24 4-watt, 254nm	95-0269-01 (115V) 95-0269-02 (230V)
Replacement Filter	38-0183-04

Accessories for EL Series UV Lamps:

- Model J-138 lamp stand for use with any EL Series Lamps allows for hands-free use
- The C-10E Cabinet (see page 20 for information) can be used with the 4-watt EL Lamps


Accessories

J-138 Lamp Stand	18-0063-01
C-10E Cabinet (for 4W EL Lamps)	95-0072-08

J-138 Stand Dim: 10.74H x 13L x 3W in.
(273 x 330 x 76mm)

J-138 Stand shown with 8-watt lamp.

EL Series Basic Model 4-watt UVS-14 includes UV and white light.



Basic Models Ordering Information

UVL-18 8-watt, 365nm/white	95-0198-01 (115V) 95-0198-02 (230V)
Replacement Filter	38-0158-04
UVM-18 8-watt, 302nm/white	95-0199-01 (115V) 95-0199-02 (230V)
Replacement Filter	38-0158-05
UVS-18 8-watt, 254nm/white	95-0200-01 (115V) 95-0200-02 (230V)
Replacement Filter	38-0158-05
UVL-16 6-watt, 365nm/white	95-0272-01 (115V) 95-0272-02 (230V)
Replacement Filter	38-0182-02
UVM-16 6-watt, 302nm/white	95-0273-01 (115V) 95-0273-02 (230V)
Replacement Filter	38-0182-01
UVS-16 6-watt, 254nm/white	95-0274-01 (115V) 95-0274-02 (230V)
Replacement Filter	38-0182-01
UVL-14 4-watt, 365nm/white	95-0264-01 (115V) 95-0264-02 (230V)
Replacement Filter	38-0183-02
UVS-14 4-watt, 254nm/white	95-0266-01 (115V) 95-0266-02 (230V)
Replacement Filter	38-0183-01

Replacement Tubes

Tube, 8-watt, White	34-0056-01
Tube, 8-watt, 365nm	34-0006-01
Tube, 8-watt, 302nm	34-0042-01
Tube, 8-watt, 254nm	34-0007-01
Tube, 6-watt, White	34-0063-01
Tube, 6-watt, 365nm	34-0034-01
Tube, 6-watt, 302nm	34-0044-01
Tube, 6-watt, 254nm	34-0013-01
Tube, 4-watt, White	34-0066-01
Tube, 4-watt, 365nm	34-0005-01
Tube, 4-watt, 254nm	34-0003-01

