

Beckman Rotor Equivalents				
	FiberLite轉子	機型	原廠轉子	最高轉速
高轉速	F9B-4x1000y	J2 Series J2-21M	JLA-9.100	8000 8000
	F10BA-6x500y	J2 Series Avanti J Series (except the Avanti J-	JLA-10.500,JA-10	10000 10000
	F10B-6x500y	J2 Series J2-21M	JA-10	10000 10000
	F14BA-6x250y*	Avanti J Series (except the Avanti J- Avanti J-E	JA-16.250	14000 14000
	F14B-6x250y	J2 Series J2-21M	JA-14	14000 14000
	F13B-14x50cy	J2 Series J2-21M	JA-12	13000 13000
	F14BA-14x50y*	Avanti J Series Avanti J-E	JA-12	14000 14000
	F18B-12x50	J2 Series J2-21M	JA-17	17000 17000
	F18BA-12x50	Avanti J Series (except the Avanti J-	JA-17	17000
	F21B-8x50y	J2 Series J2-21M	JA-20	21000 21000
超高轉速	F40L-8x100	Floor model L-Series	Type 45Ti	40000
	F50L-8x39	Floor model L-Series	Type 42.1Ti Type 55.2Ti Type 50.2Ti SW28	50000
	F65L-6x13.5	Floor model L-Series	Type 90Ti	65000
	F50L-24x1.5	Floor model L-Series	TLA-55	50000
低轉速	F8B-6x1000y	J6 J6M	N/A	6000 6000
	H3-LV	J6 J6M	N/A	3200 3200
桌上型	F14B-6x250LE	Allegra 6 Series GS-6 Series GP Series	N/A	6000 6000 6000
	F13B-14x50cy	Allegra 6 Series GS-6 Series GP Series	N/A	6000 6000 6000
	H3-LV	Allegra 6 Series GS-6 Series GP Series	G.H.3.8	3200 3200 3200

Hitachi Rotor Equivalents				
	FiberLite轉子	機型	原廠轉子	最高轉速
高轉速	F9Hi-4x1000y	CR21/22 Series	R9AF	9000
	F9Hi-6x500	CR21/22 Series	R10A R10A5	9000
	F13Hi-6x250	CR21/22 Series	R13AF3	12500
	F13Hi-14x50cy	CR21/22 Series	R15A	12500
低轉速	F7Hi-6x1000y	CR7 Series	R7AF2 R7AF3	7000

Eppendorf Rotor Equivalents				
	FiberLite轉子	機型	原廠轉子	最高轉速
桌上型	F15EP-6x100	5810 5804	F-34-6-38	12000
	F15EP-8x50c	5810 5804	F-34-6-38 (with adapters)	12000
	F21EP-48x1.5/2.0	5810 5804	F45-30-11 FA45-30-11	14000
	F21EP-30x1.5/2.0	5810 5804	F45-30-11 FA45-30-11	14000
	H4-MTP	5810 5804	A-4-62-MTP	4500
	F21EP-30x1.5/2.0	5417	F45-30-11 FA45-30-11	14000
	F21EP-24x1.5/2.0	5417	24-piace rotor	14500

Thermo Scientific Rotor Equivalents				
	FiberLite轉子	機型	原廠轉子	最高轉速
高轉速	F8S-6x1000y	Evolution RC	None	8500
	F9S-4x1000y	RC-5C Plus	None	7000
		RC-5C		7000
		RC-26 Plus		7000
		Evolution RC		9000
		RC-6		7000
		RC-6 Plus		8000
		RC-5B Series		7000
		KR22i		9000
		kr25i / HS-23		9000
	F10S-6x500y	RC-5C Plus	SLA-3000 GS-3 AK-500.11	10000
		RC-5C		9000
		RC-26 Plus		10000
		Evolution RC		12000
		RC-6		10000
		RC-6 Plus		10000
		RC-5B Series		10000
		KR22i		10000
	F14S-6x250y	RC-5C Plus	SLA-1500 GSA AK-250.14	14000
		RC-5C		13000
		RC-26 Plus		14000
		Evolution RC		14000
		RC-6		14000
		RC-6 Plus		14000
		RC-5B Series		14000
		KR22i		14000
	F20S-6x100	RC-5C Plus	SA-800	20000
		RC-26 Plus		20000
		Evolution RC		18000
		RC-6 Plus		20000
		RC-5B Series		20000
		KR22i		18000
	F21S-8x50y	RC-5C Plus	SS-34	20000
		RC-26 Plus		20000
		Evolution RC		20500
		RC-6		20000
		RC-6 Plus		18000
		KR22i		20000
		RC-5B Series		20000
	F13S-14x50cy	RC-5C Plus	SS-34 SLA-600TC	13000
		RC-5C		9000
		RC-26 Plus		13000
		Evolution RC		13000
		RC-6		13000
		RC-6 Plus		13000
		KR22i		13000
		RC-5B Series		13000
	F18S-12x50	RC-5C Plus	SLA-600TC SA-600	16500
		RC-5C		16500
		RC-26 Plus		16500
		Evolution RC		16500
		RC-6		16500
		RC-6 Plus		16500
		KR22i		16500
		RC-5B Series		16500
	F21S-48x1.5/2.0	RC-5C Plus	F-20/Micro	20000
		RC-26 Plus		20000
		Evolution		20000
		RC-6 Plus		20000
		KR22i		18000
		RC-5B Series		20000
超高轉速	F40L-8x100	Discovery Series	A-641	40000
		WX Series	T647.5	40000
	F50L-8x39	Discovery Series	A-841	50000
		WX Series	T-865 SureSpin 630/36	50000
	F65L-6x13.5	Discovery Series	T-865.1	65000
		WX Series	T-875 T-890	65000
	F50L-24x1.5	Discovery Series	N/A	50000
		WX Series		50000