



Hydrogen 19" Rack Plus

For Continuous Emission and Air Monitoring Applications

The same features found in our laboratory gas generators are now available in industry standard 19" rack mount format – top quality, ultra safe, supremely reliable, and economical hydrogen generation in a convenient and compact package.

HYDROGEN 19" RACK PLUS FEATURES:

- H₂ flow rates from 100 to 1000 ml/min
- H₂ purity options 99.9996% or 99.99996%
- H₂ pressure up to 11 bar (160 psig)
- Zero air flow rates 1.8 or 5.0 lpm options
- Zero air purity <0.1 ppm
- Minimal maintenance