

Mini UV Lamps



- Mini UV Lamps are 4-watt models available with longwave (365nm), shortwave (254nm) UV or multiband split-tube shortwave and longwave in one unit. Lamps operate on four AA batteries (not included).

Mini UV Lamps
UVG-4, UVSL-14P,
UVL-4

- Excellent lamps for office or field inspection
- Lamps are lightweight and small enough to fit into a pocket
- A nylon strap is attached for easy handling

Dimensions: **UVG-4, UVSL-14P** 6.25L x 1W x 2.5H in. (159 x 25 x 64mm) **UVL-4** 6.25L x 1W x 2H in. (159 x 25 x 51mm); Weight: 0.5 lb (0.23 kg)

Rechargeable UV Lamps

Rechargeable 6-watt 26P series lamps come in three dual-tube models.

- Lamps are equipped with an internal gel based, lead-acid battery which can be charged for field operation from a standard wall outlet
- Optional 12V adapter is available for operating the lamp from a car cigarette lighter
- The UVSL-26P (shown) can operate with both long and shortwave on at the same time
- Use with the C-10P Cabinet



Rechargeable UV Lamps with built-in rechargeable battery.

Dimensions: 5.5L x 3.0W x 9.8H in. (249 x 76 x 140mm); Weight: 4 lb (1.8 kg)

Rechargeable Lamps - Ordering Information

| | |
|--|--|
| UVL-26P 365nm Replacement Tube | 95-0186-01 (115V) 95-0186-02 (230V) 34-0154-02 |
| UVS-26P 254nm Replacement Tube Replacement Filter | 95-0187-01 (115V) 95-0187-02 (230V) 34-0013-01 38-0154-01 |
| UVSL-26P 254nm/365nm Replacement Tube, 365nm Replacement Tube, 254nm Replacement Filter Adapter, 12V | 95-0181-01 (115V) 95-0181-02 (230V) 34-0016-01 34-0013-01 38-0154-02 58-0127-02 |




UVL-4F Mini UV Lamp
includes a longwave UV tube
and a flashlight.

Dimensions UVL-4F:
7L x 1.3W x 2.5H in. (178 x 33 x 64mm)
Weight: 0.5 lb (0.23kg)

Ordering Information

| | |
|---|--------------------------------|
| UVG-4 254nm Replacement Tube | 95-0158-04 (4AA) 34-0003-01 |
| UVL-4 365nm BLB Replacement Tube | 95-0125-05 (4AA) 34-0010-01 |
| UVSL-14P Multiband split-tube 254nm/365nm Replacement Tube | 95-0188-02 (4AA) 34-0004-01 |
| UVL-4F 365nm BLB/Flashlight Replacement Tube | 95-0305-01 (4AA) 34-0010-01 |

Portable UV Lamps



Portable UV Lamps
lightweight for field use.

Portable 6-watt UV lamps are operated by two 6-volt lantern batteries.

- The distinctive wrap around handle comfortably balances the lamp in your hand
- A built-in flashlight is included
- Optional rechargeable 6V batteries and charger are available

Dimensions: 9.5L x 2.8W x 9.4H in. (241 x 71 x 238mm); Weight: 4 lb (1.8 kg)

Portable Lamp - Ordering Information

| | |
|---|--|
| ML-49 365nm BLB Replacement Tube | 95-0011-01 (12V) 34-0016-01 |
| UVG-47 254nm Replacement Tube Replacement Filter | 95-0009-02 (12V) 34-0013-01 98-0013-01 |
| UVGL-48 split-tube 254nm/365nm Replacement Tube Replacement Filter Bulb, White Battery, 6V Battery, rechargeable 6V (115V) Charger (115V) | 95-0010-02 (12V) 34-0015-01 98-0010-04 34-0014-01 45-0005-01 45-0005-02 58-0136-01 |