



Mistral Evolution



Whisper

## Mistral Evolution and Whisper Series

### For LC/MS/MS and ELSD Applications

The Mistral Evolution and Whisper nitrogen generators have been specifically designed to meet the flow, pressure and purity requirements of all the leading LC/MS instrument manufactures. Nitrogen is produced by either carbon molecular sieve (CMS) or hollow fiber membrane technology, depending on the exact user requirements. Generators are also available with integral compressors, providing a complete "plug and play" solution.

#### **MISTRAL EVOLUTION AND WHISPER SERIES FEATURES:**

- Flow rates from 10 to 120 lpm
- LC/MS grade hydrocarbon free nitrogen
- Pressure up to 8 bar (116 psig)
- Oil free, low noise compressor options
- Specialist hybrid air and nitrogen options for Agilent Microwave Plasma and Sciex LC/MS/MS