



## NM Plus Series

**Ideal for all GC & GC/MS carrier gas and detector gas applications**

The NM Plus Series enhances the high performance of the PG Plus with the addition of a unique no maintenance purification system. This results in ultra high purity carrier grade gas. With its highly sophisticated software control and safety alarm capability it presents the GC user with the opportunity to reap the benefits of hydrogen carrier gas while overcoming the safety concerns.

### NM PLUS SERIES FEATURES:

- Flow rate from 100 to 1000 ml/min
- Purity 99.99996%
- Pressure up to 11 bar (160 psig)
- No maintenance
- USB for remote monitoring and servicing
- Unique 100% titanium electrolytic cell
- 2 year standard warranty