



Nitrogen Generators

For Gas Chromatography and Analytical Instrument Applications

VICI DBS nitrogen generators are engineered to transform standard compressed air into a safe regulated supply of Ultra High Purity nitrogen. Utilizing a superior quality carbon molecular sieve (CMS) and innovative column switching, the nitrogen generators require minimal maintenance and operator attention.

NITROGEN GENERATOR FEATURES:

- Flow rates from 500 to 4,000 ml/min
- Purity up to 99.999%
- Pressure up to 5 bar (75psig)
- Space saving design
- Optional hydrocarbon removal for GC applications
- Ideal for GC carrier gas including ECD, make up gas, TOC, DSC, and ELSD applications