

冷凍凝結槽

Cold Traps

UT 系列



選配附件




UT-1000 用上蓋 B



UT-2000 用上蓋 B2



酸、有機溶劑用 Set D (500 ml) 低沸點、有機溶劑用 Set D (500 ml)




水系 Trap


有機溶劑系 Trap



UT-3000A・4000A 用



UT-3000L・4000L 用



UT-3000A・4000A 用
UT-3000L・4000L 用



酸、有機溶劑用 Set C (500 mL) UT-3000A・4000A 用



大容量酸、有機溶劑用 Set C (1800 mL) UT-3000L・4000L 用




低沸點有機溶劑用 Set C (500 mL) UT-3000A・4000A 用



搭配真空烘箱



搭配離心濃縮機



搭配真空旋轉濃縮機

型號	UT-1000	UT-2000	UT-3000A	UT-3000L	UT4000A	UT-4000L
除溼容量 (水)	max. 0.5kg	max. 0.4kg	max. 1kg	max. 4kg	max. 1kg	max. 4kg
最低到達溫度	-50 °C	-100 °C	-50 °C		-80 °C	
冷凍機功率	300W	300W+300W	450W	675W	450W+450W	450W+450W
內槽尺寸 (mm) ; 容量	內徑 90 x 225 高, 約 1.4 L		直徑 153 x 225 高, 約 4 L	直徑 200 x 395 高, 約 9.4 L	直徑 153 x 255 高, 約 4 L	直徑 200 x 395 高, 約 9.4 L
外部尺寸 (mm)	373W x 370D x 317H	461W x 510D x 473H	335W x 362D x 838H	372W x 382D x 835H	362W x 482D x 838H	377W x 532D x 835H

上蓋及 Glass Condenser 為選配之附件