

# 實驗室型凍結乾燥機

Freeze Dryers

**FD-1000**

**FDU-1200 · 2200**



-45 °C  
500 mL /batch

FD-1000



-45 °C  
1 L /batch

FDU-1200

-80 °C  
1 L /batch



FDU-2200



可接 filter 防污染




FD-1000 操作面板



FDU-1200 · 2200 操作面板

## ■ 組合例



型號	FD-1000	FDU-1200	FDU-2200
冷卻溫度	-45 °C	-45 °C	-80 °C
乾燥能力	500 ml/batch	1 L/batch	1 L/batch
冷凍槽尺寸 (mm)	φ 83.1 x 277H	φ 200 x 175H	φ 153 x 225H
冷凍機功率	250W	350W	350W x 2
外部尺寸 (mm)	195W x 440D x 790H	400W x 420D x 440H	450W x 550D x 935H

\*需另配合真空幫浦使用