


Grinding/Polishing

Automatic grinding head/tip can be changed at both sides for use.



P20FS-1-A3; P20FS-1-A5
P25FS-1-A3; P25FS-1-A5



P20FRS-A3; P20FRS-A5
P25FRS-A3; P25FRS-A5

MODEL	P20FS-1-A3/A5		P25FS-1-A3/A5		P20FRS-A3/A5		P25FRS-A3/A5	
No. of discs	1				2			
Grinding Method	Automatic							
Disc dia	φ 200		φ 250		φ 200		φ 250	
Rotational speed	0 ~ 600rpm				0 ~ 600rpm			
Motor power	DC180V or DC90V 1φ 180w				DC180V or DC90V 1φ 250w			
Cooling unit	1 piece				2 piece			
Machine type	Table top				Table top			
No. of specimens	A3/Max. 3 Pcs; A5/Max. 5 Pcs;				A3/Max. 3 Pcs; A5/Max. 5 Pcs;			
Workpiece	32mm (special size, optional φ20, φ25, φ30)							
Force per specimen	0 ~ 5kg adjustable							
Rotational speed	120 rpm fixed							
Timer control	Liquid crystal display, the grinding time setting is available.							
Dimensions	740(W) x 440(D) x 540(H)mm				740(W) x 700(D) x 540(H)mm			
Power	AC 220V 1φ or AC 110V 1φ				AC 220V 1φ or AC 110V 1φ			
Weight	52kg		54kg		70kg		75kg	
Standard accessories	Aluminum Oxide Powder		100g		Aluminum Oxide Powder		100g	
	Drain Pipe	1.5M	Anti-splash Cover	1pc	Drain Pipe	1.5M	Anti-splash Cover	2pc
	Small Pipe Clip	2pcs	Water Sand Paper	1pc	Small Pipe Clip	2pcs	Water Sand Paper	2pc
	Abrasive Cloth	1pc	Big Pipe Clip	2pcs	Abrasive Cloth	2pc	Big Pipe Clip	2pcs
	Water Inlet Pipe	1.5M	O Ring	1pc	Water Inlet Pipe	1.5M	O Ring	1pc
Accessories for optional purchase	W1 Automatic drip bottle function							

Agent