

# 小型多檢體有機合成裝置

Personal Organic Synthesizers

## Zodiac

### CCX-1000 · 1200


強力還流攪拌  
CCX-1000



-78 °C 冷卻實驗  
CCX-1200



型式	CCX-1000	CCX-1010	CCX-1200	CCX-1210
試驗管數	1 ~ 12 支	1 ~ 12 支	1 ~ 12 支	1 ~ 12 支
合成量	3 ~ 15 mL	3 ~ 15 mL	3 ~ 15 mL	3 ~ 15 mL
使用溫度	室溫+10 ~ 160 °C	室溫+10 ~ 160 °C	-78 °C ~ 室溫	-78 °C ~ 室溫
控溫方式	冷卻水	冷卻水	乾冰+丙酮	乾冰+丙酮
還流	√	√	×	×
標準配件	試驗管 x 12, 攪拌子 x 12, Septum x 12, 溫度 Sensor x 1			



強力攪拌




密閉簡單



試葯添加



-78 °C 冷卻實驗



冷卻循環 (CCX-1000)



個別取用



移動方便



個別氣體置換